

# Roche LightCycler 480 Real-time PCR平台

## 原理暨應用介紹



## LIGHTCYCLER® 480 REAL-TIME PCR SYSTEM

The LightCycler® 480 System is a proven high-performance, medium-to high-throughput PCR platform that provides various methods for gene detection, gene expression analysis, genetic variation analysis, and array data validation.

*Sophia Lin*

又鑫生物科技有限公司

Real Time PCR Basic Training

Troubleshooting cases sharing

System Operation Procedures

LC 480 QC Report

Q&A

Real Time PCR Basic Training



Troubleshooting cases sharing

System Operation Procedures

LC 480 QC Report


Q&A

# Polymerase Chain Reaction (PCR)




# Real-Time Polymerase Chain Reaction (real-time PCR)

a SYBR Green assay



b TaqMan assay








**Non-specific binding dye**



**Specific labeled probes**

SYBR Green I







Hybridization Probe

Hydrolysis Probe (TaqMan)






# Analysis Module **T<sub>m</sub> Calling/Melting Curve**









## *Fluorescence Resonance Energy Transfer (FRET)*



# LightCycler® Assay Formats

## TaqMan Probe




1. 多色Multiplex
2. 不需melting curve




# LightCycler® Assay Formats

## *Dye Emission Spectra*



# LightCycler® Assay Formats

## *Dye Emission Spectra and Crosstalk*



# Tricolor Hydrolysis Probe - Example

## *Application with Spectral Crosstalk*

FAM


—

Red 610


—

Cy 5 (PBGD + G6PDH + CyP2 C9)


Crosstalk of Cy 5 into channel 610 is reliably compensated




533



610



670 no CC




670 with CC

# LightCycler® Assay Formats


## *Hybridization Probes*

HybPr<sub>be</sub>


**A** Denaturation



**B** Annealing



**C** Elongation



**D** End of Elongation




No primer-dimer or side-products' disturbance ! – High specificity

Don't need to run the Melting curve – High efficiency,


# There is a correlation between Cp and concentration

*The **higher concentration** of target nucleic acid in the starting material, the sooner a significant increase in fluorescent signal will be observed, yielding a **lower Cycle no.***



# PCR Quantification

## *Theoretical and Practical Aspects*



log-phase-PCR



end-point-PCR

N: number of amplified molecules





$N_0$ : initial number of molecules

n: number of amplification cycles


E: amplification efficiency



# How to calculate the Ct/Cp/Cq value






# Concentrations and Crossing Points






## Calculation of crossing points (Cp)

- Optional with **Fit Point** method (user-influenced)
- Standard method (automatic)  
**2nd Derivative Maximum Method**  
二次微分最大値

# Fit Point method (user-influenced)




# Fit Point method (user-influenced)




# Crossing Point Calculation

## *2<sup>nd</sup> Derivative Maximum Method*




# There is a correlation between $C_p$ and concentration


The **higher concentration** of target nucleic acid in the starting material, the sooner a significant increase in fluorescent signal will be observed, yielding a **lower Cycle no.**



# Absolute Quantification with External Standards: Principle



# Absolute Quantification with External Standards: Example





# Relative Quantification *Scheme*

**Wanted:**


Target gene expression level

**Required:**

Target gene concentration

Reference gene concentration

Calibrator sample (optional)



# Relative Quantification - Without Efficiency correction

$$\text{Relative Ratio} = \frac{\text{concentration of target}}{\text{concentration of reference}}$$

| Type   | gene      | Cp |
|--------|-----------|----|
| Sample | Target    | 22 |
|        | Reference | 21 |

$$\begin{aligned} \frac{22}{21} &\rightarrow \frac{2^{22}}{2^{21}} \\ &= 2^{22-21} = 2^1 \rightarrow = 2^{-1} = 1/2 \\ &= 2^{\Delta Ct} \rightarrow = 2^{-\Delta Ct} \end{aligned}$$

| Type                    | gene      | Cp |
|-------------------------|-----------|----|
| Calibrator<br>(control) | Target    | 25 |
|                         | Reference | 20 |

$$\begin{aligned} &= 2^{25-20} = 2^5 \rightarrow = 2^{-5} \\ &= 2^{\Delta Ct} \rightarrow = 2^{-\Delta Ct} = 1/32 \end{aligned}$$

# Relative Quantification - Without Efficiency correction

$$\text{Calibrator Normalized Ratio} = \frac{\frac{\text{concentration of target (sample)}}{\text{concentration of reference (sample)}}}{\frac{\text{concentration of target (calibrator)}}{\text{concentration of reference (calibrator)}}}$$

| Type   | gene      | Cp |
|--------|-----------|----|
| Sample | Target    | 22 |
|        | Reference | 21 |

$$= 2^{22-21} = 2^1 \rightarrow = 2^{-1} = 1/2$$

$$= 2^{\Delta Ct} \rightarrow = 2^{-\Delta Ct}$$

| Type                    | gene      | Cp |
|-------------------------|-----------|----|
| Calibrator<br>(control) | Target    | 25 |
|                         | Reference | 20 |

---


$$= 2^{25-20} = 2^5 \rightarrow = 2^{-5} = 1/32$$

$$= 2^{\Delta Ct} \rightarrow = 2^{-\Delta Ct}$$

$$= 2^{-1-(-5)} = 2^4 = 16$$

$$= 2^{-\Delta\Delta Ct}$$

# Example for Relative Quantification

| Type       | gene      | Cp |
|------------|-----------|----|
| Calibrator | Target    | 25 |
|            | Reference | 20 |
| Sample     | Target    | 22 |
|            | Reference | 21 |

## I. Without efficiency correction

$$\begin{aligned}
 \text{Relative amount} &= 2^{-\Delta\Delta CT} \\
 &= 2^{-[(22-21)-(25-20)]} = 2^{-[1-5]} = 2^4 = 16
 \end{aligned}$$




## II. With efficiency correction

$$\begin{aligned}
 &\xrightarrow{\quad} \begin{aligned} E_T &= 1.6 \\ E_R &= 1.8 \end{aligned}
 \end{aligned}$$

$$\begin{aligned}
 \text{Relative amount} &= E_T^{CpT(C) - CpT(S)} \times E_R^{CpR(S) - CpR(C)} \\
 &= 1.6^{(25-22)} \times 1.8^{(21-20)} = 1.6^3 \times 1.8^1 = 7.37
 \end{aligned}$$





# Endpoint Genotyping (allelic discrimination)

## Principle





# Endpoint Genotyping (allelic discrimination)

## Principle



# Endpoint Genotyping (allelic discrimination)




## Principle



# Endpoint Genotyping Data

## *Influence of DNA Concentration*




1 and 20 ng DNA






# Melting Curve Genotyping

## *Principle*



# Melting Curve Genotyping

## *Principle*




Low

Medium


High

Temperature



# Melting Curve Based Genotyping



## *Genotyping of Single Point Mutation*



# Endpoint & Melting Genotyping Data



## *Influence of DNA Concentration*

1 and 20 ng DNA



# Melting Genotyping Data

## *SimpleProbe*



# NAT 檢測流程



## Real-Time PCR Application

- Absolute Quantification
- Relative Gene Expression Analysis
- Diagnostic / Biomarker Analysis
- Genotyping / Mutation Detection
- Environment Microbial Detection
- Food Testing
- Protein Shift Assay

Real Time PCR Basic Training

**Troubleshooting cases sharing**

System Operation Procedures

LC 480 QC Report





Q&A

## Troubleshooting cases sharing:

- **Shoulder of Melting Curve**
- **Primer Dimmer 1**
- **Primer Dimmer 2**
- **gDNA contamination**
- **Effect of master mix**




# Case 1: Shoulder of Melting Curve




## Case 2: Primer Dimmer

### *No Template Control*



Positive Sample Melting curve




#### 1. Primer dimer :

Primer conc. ↓

Primer design

Amplicon length : 100bps↑




#### 2. Contamination


## Case 3: Primer Dimmer

*Not all primer dimers are a problem for an Assay*

NTC



Sample Results



# How to Design Primers for QPCR




- PCR product/Amplicon size : 50-210 base pairs.
- Primer length : 19-23 nucleotides.
- GC content : 35-65%.
- Melting temperature ( $T_m$ ) : 60-68 °C. The annealing temperature for the assay is 5 °C lesser than the  $T_m$  of the primers.
- Exon-exon junction : When amplifying cDNA by QPCR, the primers should span exon-exon junction to avoid the amplification of contaminating DNA.
- Repeats and runs : Dinucleotide repeats ( TCTCTCTCTC) and repeated nucleotides (eg. TAAAAAAGC) should be avoided.
- 3' Complementarity : The complementary regions of the 3' ends of forward and reverse primers should be avoided to prevent the formation of primer-dimers.
- 3' Stability : G or C residues should be included at the 3' end of the primer to increase the stability of the annealing.
- GC clamp : One or two GC clamps at the 5' end of the primer increases the specificity of the annealing.
- Specificity : The specificity of the primers should be checked by BLAST
- SNPs : Primers should not contain any known SNP (single nucleotide polymorphism) variations




# Case 4: gDNA contamination

## Cross Intron Assay Design



Sample Results



No Reverse Transcription



# Case 5: Effect of Master Mix



# Tips for successful Real-Time PCR

- Avoid of contamination.
- High quality RNA.
- Good primer and probe design.
- Decrease popettes error.
- Stable master mix and mix well all reagent.
- Control: PC, NC, NTC, NRT.
- Avoid primer dimer: melting curve.

Real Time PCR Basic Training

Troubleshooting cases sharing

**System Operation Procedures**

LC 480 QC Report

Q&A




# LightCycler<sup>®</sup> 480 Real-Time PCR System



# Data Capture (Other Brand)

Fig. 12



**Report signal : SYBR Green 1**

**Reference signal : ROX**


• 校正

**(MeltDoctor™ HRM Calibration Plates)**

# Data Capture (Other Brand)




- 非同步
- 時間長



- 校正  
(Melt Calibration Kit)

# LightCycler® 480 System

## *Optical unit and lightpath*



- **LED lamp**  
high intensity  
broad dynamic range  
lifetime  
for **LED**: approx. 10,000 hrs
- **CCD camera**
- **Five excitation filters**
- **Six detection filters**
- **Optimized arrangement of optical components**
- **Homogeneous excitation and fluorescence detection**

Filter Combination Selection

|            |     | Emission                 |                                     |                                     |                          |                          |                                     |
|------------|-----|--------------------------|-------------------------------------|-------------------------------------|--------------------------|--------------------------|-------------------------------------|
|            |     | 488                      | 510                                 | 580                                 | 610                      | 640                      | 660                                 |
| Excitation | 440 | <input type="checkbox"/> | <input type="checkbox"/>            | <input type="checkbox"/>            | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>            |
|            | 465 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/>            | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>            |
|            | 498 | <input type="checkbox"/> | <input type="checkbox"/>            | <input type="checkbox"/>            | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>            |
|            | 533 | <input type="checkbox"/> | <input type="checkbox"/>            | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>            |
|            | 618 | <input type="checkbox"/> | <input type="checkbox"/>            | <input type="checkbox"/>            | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> |
|            |     |                          |                                     |                                     |                          |                          |                                     |

Clear

---

Selected Filter Combination List

| Excitation Filter | Emission Filter | Name            | Melt Factor | Quant Factor | Max Integration Time (Sec) |
|-------------------|-----------------|-----------------|-------------|--------------|----------------------------|
| 465               | 510             | FAM             | 1           | 10           | 2                          |
| 533               | 580             | VIC / HEX / Y1  |             | 10           | 2                          |
| 618               | 660             | Cy 5 / Cy 5.5 1 |             | 10           | 2                          |

# System Start-Up



| 左 警示燈  | 右 警示燈  | 機器狀態                 |
|--------|--------|----------------------|
| 橘 *閃爍* | 橘 *閃爍* | 正在初始化                |
| 綠      | 橘      | 機器啟動完成，96/384 孔盤還未放入 |
| 綠      | 橘 *閃爍* | 96/384 孔盤正在放入中       |
| 綠      | 綠      | 機器啟動完成，96/384 孔盤已經放入 |
| 綠 *閃爍* | 綠 *閃爍* | 實驗進行中                |

# System Start-Up

- ▶ To load the prepared multiwell plate into the LightCycler<sup>®</sup> 480 Instrument, press the push button on the front of the instrument (located next to the instrument status LEDs):



Place the multiwell plate into the loading frame of the loader with the flat edge pointing towards the instrument. (The short plate edge with beveled corners points away from the instrument.)



Press the plate loading push button again to retract the loader with the inserted multiwell plate into the instrument. You are now ready to start the run.





# LightCycler® 480 System

## *Disposables*




96-well plate for  
10–100  $\mu$ l



LightCycler® 8-Tube Strips (white and clear)

# Open the software and Login to the data base



Instrument: No active instrument Database:

Window:  User:




For PC login  
ID: **ICOB**  
PW: **non**

Login

User name: \*  admin

Password: \*  Lightcycler480

Log on to: \*


 Options >>  

56



# LightCycler® 480 Software Version 1.5

## Front screen: New experiment from template



The screenshot displays the LightCycler® 480 Software release 1.5.0.0804 interface. The main window shows the 'Run Templates' list, which is a table of templates. The 'Gene Scanning 96-II' template is selected. On the right side, the 'Experiment Creation' panel is visible, with the 'New Experiment from Template' button highlighted by a red rectangle.

**Run Templates**

| Name  | Path                               | Creation Date           |
|---|------------------------------------|-------------------------|
| Dual Color Hydrolysis Probe - UPL Probe 96-II | /Roche/Templates/Run Templates/Sys | 07/20/2007 16:33:50.612 |
| Endpoint Genotyping (PCR Read) 96-II          | /Roche/Templates/Run Templates/Sys | 07/20/2007 17:08:36.713 |
| Endpoint Genotyping (Pre-Post Read) 96-II     | /Roche/Templates/Run Templates/Sys | 07/20/2007 17:13:51.398 |
| Gene Scanning 96-II                           | /Roche/Templates/Run Templates/Sys | 07/20/2007 16:44:38.871 |
| HybProbe 96-II                                | /Roche/Templates/Run Templates/Sys | 07/18/2007 15:44:17.911 |
| Mono Color Hydrolysis Probe - UPL Probe 96-II | /Roche/Templates/Run Templates/Sys | 07/20/2007 16:32:53.570 |
| SimpleProbe 96-II                             | /Roche/Templates/Run Templates/Sys | 07/20/2007 15:53:27.964 |
| SYBR Green I 96-II                            | /Roche/Templates/Run Templates/Sys | 07/20/2007 15:54:55.591 |
| 20160310                                      | /System Admin/Templates/Run Templa | 03/10/2016 16:02:07.684 |
| Test Run Protocol                             | /System Admin/Templates/Run Templa | 03/08/2016 11:40:51.917 |

**Experiment Creation**

Plate Type

- ☒ White Plates
- ☐ Clear Plates

New Experiment

New Experiment from Macro

**New Experiment from Template**

Open Existing Object

**Absolute Quantification**

# Start Run



Experiment

Subset Editor

Sample Editor

Analysis

Report

Run Protocol

Data

Run Notes

Setup

Detection Format SYBR Green I / HI Customize
Block Size 96
Plate ID 
Reaction Volume 20

Color Comp ID 
Lot No 
Test ID

Programs

| Program Name   | Cycles | Analysis Mode  |
|----------------|--------|----------------|
| pre-incubation | 1      | None           |
| amplification  | 45     | Quantification |
| melting curve  | 1      | Melting Curves |
| cooling        | 1      | None           |

Temperature Targets

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |
|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| 95          | None             | 00:05:00        | 4.4              |                       | 0               | 0              | 0                   |

Overview

Apply Template

Save As Template

End Program

+ 10 Cycles

Start Run

# Start Run



Experiment

Subset Editor

Sample Editor

Analysis

Report

| Run Protocol  |                  | Data  |                  | Run Notes                                   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
|---|------------------|---|------------------|---|-----------------|----------------|---------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|----------------|---------------|----------|----------------|---------|---|------|---|----|------|----------|-----|--|---|---|---|----|--------|----------|-----|--|---|---|---|
| <div> <div>Setup</div> <div> <div>Detection Format</div> <div>SYBR Green I /</div> <div>Customize</div> </div> <div> <div>Block Size</div> <div>96</div> </div> <div> <div>Plate ID</div> <div></div> </div> <div> <div>Reaction Volume</div> <div>20</div> </div> </div>   |                  |   |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| <div> <div>Color Comp ID</div> <div></div> </div>   |                  | <div> <div>Lot No</div> <div></div> </div>                                      |                  | <div> <div>Test ID</div> <div></div> </div> |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| <div> <div>↑</div> <div>Programs</div> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th>Program Name</th> <th>Cycles</th> <th>Analysis Mode</th> </tr> </thead> <tbody> <tr> <td>pre-incubation</td> <td>1</td> <td>None</td> </tr> <tr> <td>amplification</td> <td>45</td> <td>Quantification</td> </tr> <tr> <td>melting curve</td> <td>1</td> <td>Melting Curves</td> </tr> <tr> <td>cooling</td> <td>1</td> <td>None</td> </tr> </tbody> </table> <div>↓</div> </div>   |                  |   |                  |   |                 | Program Name   | Cycles              | Analysis Mode   | pre-incubation   | 1                     | None            | amplification  | 45                  | Quantification | melting curve | 1        | Melting Curves | cooling | 1 | None |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| Program Name  | Cycles           | Analysis Mode   |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| pre-incubation  | 1                | None  |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| amplification   | 45               | Quantification  |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| melting curve   | 1                | Melting Curves  |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| cooling   | 1                | None  |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| <div> <div>↑</div> <div>amplification Temperature Targets</div> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th>Target (°C)</th> <th>Acquisition Mode</th> <th>Hold (hh:mm:ss)</th> <th>Ramp Rate (°C/s)</th> <th>Acquisitions (per °C)</th> <th>Sec Target (°C)</th> <th>Step Size (°C)</th> <th>Step Delay (cycles)</th> </tr> </thead> <tbody> <tr> <td>95</td> <td>None</td> <td>00:00:10</td> <td>4.4</td> <td></td> <td>0</td> <td>0</td> <td>0</td> </tr> <tr> <td>60</td> <td>None</td> <td>00:00:10</td> <td>2.2</td> <td></td> <td>0</td> <td>0</td> <td>0</td> </tr> <tr> <td>72</td> <td>Single</td> <td>00:00:10</td> <td>4.4</td> <td></td> <td>0</td> <td>0</td> <td>0</td> </tr> </tbody> </table> <div>↓</div> </div> |                  |   |                  |   |                 | Target (°C)    | Acquisition Mode    | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) | 95             | None          | 00:00:10 | 4.4            |         | 0 | 0    | 0 | 60 | None | 00:00:10 | 2.2 |  | 0 | 0 | 0 | 72 | Single | 00:00:10 | 4.4 |  | 0 | 0 | 0 |
| Target (°C)   | Acquisition Mode | Hold (hh:mm:ss)   | Ramp Rate (°C/s) | Acquisitions (per °C)                       | Sec Target (°C) | Step Size (°C) | Step Delay (cycles) |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| 95  | None             | 00:00:10  | 4.4              |   | 0               | 0              | 0                   |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| 60  | None             | 00:00:10  | 2.2              |   | 0               | 0              | 0                   |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| 72  | Single           | 00:00:10  | 4.4              |   | 0               | 0              | 0                   |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| <div> <div>Overview</div> </div>  |                  |   |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |
| <div> <div>Apply Template</div> <div>↓</div> </div>   |                  | <div> <div>End Program</div> <div>+ 10 Cycles</div> <div>Start Run</div> </div> |                  |   |                 |                |                     |                 |                  |                       |                 |                |                     |                |               |          |                |         |   |      |   |    |      |          |     |  |   |   |   |    |        |          |     |  |   |   |   |

# Start Run

**Experiment**

**Subst Editor**

**Sample Editor**

**Analysis**

**Report**

**Run Protocol**

**Setup**

Detection Format: SYBR Green I / HI **Customize** Block Size: 96 Plate ID: [ ]

Color Comp ID: [ ] Lot No: [ ] Test I: [ ]

**Programs**

- Program Name
- pre-incubation
- amplification
- melting curve
- cooling

**Temperature Targets**

| Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) |
|-------------|------------------|-----------------|------------------|-----------------------|
| 95          | None             | 00:05:00        | 4.4              |                       |

**Overview**

0:00:00 0:05:34 0:13:59 0:22:24 0:30:53 0:39:18

Estimated Time (h:mm:ss)

**Save Template**

- 123
- System Admin
  - Experiments
  - Macros
  - Preferences
  - Special Data
  - Templates
    - Analysis Templates
    - Report Templates
    - Run Templates**
- 0306
- 0409
- 111
- 123
- 20131001
- 20131024
- 20131108
- New Experiment Run Protocol
- Sample Templates


Name: Give me a name

**Apply Template**

**Save As Template**

**End Program** **+ 10 Cycles** **Start Run**

# Start your experiment by your template



## **LightCycler480 Software:**

- 1. New experiment from template**
- 2. Subset & sample set up**
- 3. Data analysis & report**

# Subset



Experiment

Subset Editor

Sample Editor

Analysis

Report

Step 2

Subsets

| ID | Name        | Analysis                            | Report                              |
|----|-------------|-------------------------------------|-------------------------------------|
| 1  | All Samples | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |

All Samples settings

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |

Step 1

+

-

Copy

Rename

Step 3

Apply

Clear

Cancel

## Sample Editor – AQ example

**Experiment**

**Subset Editor**

**Sample Editor**

**Analysis**

**Report**

**Step 1: Select Workflow**

☒ Abs Quant
 ☐ Rel Quant
 ☐ Scanning
 ☐ Color Comp  
☐ Tm
 ☐ Melt Geno
 ☐ Endpt Geno

**Step 2: Select Samples**

Subset: All Samples

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| B | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| C | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| D | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| E | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| F | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| G | U | U | U | U | U | U | U | U | U | U  | U  | U  |

**Step 3: Edit Abs Quant Properties**

Sample Name

Sample Type

☒ Unknown
 ☐ Negative Control  
☐ Positive Control/Calibrator  
☐ Standard Concentration

Auto Std Curve

Make Replicates

Auto Replicate

Clear Replicates

## RQ example

**Step 1: Select Workflow**

☐ Abs Quant
 ☒ Rel Quant
 ☐ Scanning
 ☐ Color Comp  
☐ Tm
 ☐ Melt Geno
 ☐ Endpt Geno

**Step 2: Select Samples**

Subset: All Samples

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A | U | U | U | U | U | U | U | U | U | U  | U  | U  |
| B | U | U | U | U | U | U | U | U | U | U  | U  | U  |

**Step 3: Edit Rel Quant Properties**

Sample Name

Sample Type

☒ Unknown
 ☐ Negative Control  
☐ Positive Control/Calibrator  
☐ Standard Concentration

Auto Std Curve

Gene target

Target name

Eff 2.00

☐ Target
 ☐ Reference
 ☒ Unassigned


Make Replicates

Auto Replicate

Clear Replicates




# Analysis



## 分析模式簡要說明：

| 分析模式                             | 說明                                       |
|----------------------------------|--|
| Abs Quant / 2nd Derivative Max   | 絕對定量 (2 <sup>nd</sup> Max)               |
| Abs Quant / Fit Points           | 絕對定量 (Fit Point)                         |
| Advanced Relative Quantification | 進階相對定量(2 <sup>nd</sup> Max or Fit Point) |
| Basic Relative Quantification    | 基礎相對定量 (Fit Point only)                  |
| Endpoint Genotyping              | 基因分型(適用於TaqMan probe)                    |
| Melt Curve Genotyping            | 基因分型(適用於Hyprobe)                         |
| Tm Calling                       | Tm分析                                     |

# Analysis



Experiment

Subset Editor

Sample Editor

Analysis

Report

Create New Analysis

- Abs Quant/2nd Derivative Max
- Abs Quant/Fit Points
- Advanced Relative Quantification
- Basic Relative Quantification
- Color Compensation
- Endpoint Genotyping
- Gene Scanning
- Melt Curve Genotyping
- Tm Calling

Create new analysis

Analysis Type \* Abs Quant/2nd Derivative Max

Subset \* All Samples

Program \* Amplification

Name \* Abs Quant/2nd Derivative Max for All

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| I |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| J |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| K |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| L |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| M |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| N |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| O |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| P |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |

☒ ☐

# AQ Analysis

Experiment

Subset Editor

Sample Editor

Analysis

Report

Subset: Quantification

Abs Quant results

☒ Positive
 ☒ Negative
 ☒ Standard

| Samples                             |       |                       | Results |               |          |       |
|-------------------------------------|-------|-----------------------|---------|---------------|----------|-------|
| Include                             | Color | Pos Name              | Cp      | Concentration | Stand... | Statu |
| <input checked="" type="checkbox"/> |       | C8 no template contrc |         |               |          |       |
| <input checked="" type="checkbox"/> |       | C9 negative control   |         |               |          |       |
| <input checked="" type="checkbox"/> |       | C10 Standard 1E1      |         |               | 1.00E1   |       |
| <input checked="" type="checkbox"/> |       | C11 Standard 1E2      |         |               | 1.00E2   |       |
| <input checked="" type="checkbox"/> |       | C12 Standard 1E3      |         |               | 1.00E3   |       |
| <input checked="" type="checkbox"/> |       | C13 Standard 1E4      |         |               | 1.00E4   |       |
| <input checked="" type="checkbox"/> |       | C14 Standard 1E5      |         |               | 1.00E5   |       |

Apply Template

Notes

Calculate

Subset: Quantification

Abs Quant results

☒ Positive
 ☒ Negative
 ☒ Standard

| Samples                             |       |                       | Results |               |          |       |
|-------------------------------------|-------|-----------------------|---------|---------------|----------|-------|
| Include                             | Color | Pos Name              | Cp      | Concentration | Stand... | Statu |
| <input checked="" type="checkbox"/> |       | C8 no template contrc |         |               |          |       |
| <input checked="" type="checkbox"/> |       | C9 negative control   | 35.39   | 2.82E1        |          |       |
| <input checked="" type="checkbox"/> |       | C10 Standard 1E1      | 37.45   | 5.49E0        | 1.00E1   |       |

Replicate Statistics

| Samples    | MeanCp | STD Cp | Mean conc | STD conc |
|------------|--------|--------|-----------|----------|
| C8, D8, E8 |        |        |           |          |
| C9, D9, E9 | 36.31  | 0.93   | 1.61E1    | 1.11E1   |

Apply Template

Notes

Calculate

# Analysis

Experiment

Subset  
EditorSample  
Editor

Analysis

Report

## Create New Analysis

Abs Quant/2nd Derivative Max

Abs Quant/Fit Points

Advanced Relative Quantification

Basic Relative Quantification

Color Compensation

Endpoint Genotyping

Gene Scanning

Melt Curve Genotyping

Tm Calling

## Create new analysis

Analysis Type

\* Abs Quant/Fit Points

Subset

\* All Samples

Program

\* amplification

Name

\* Abs Quant/Fit Points for All Samples

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |



# AQ Analysis


Experiment

Subset Editor


Sample Editor

Analysis

Report



# AQ Analysis



# AQ Analysis


Experiment

Subset Editor

Sample Editor

Analysis

Report



# AQ Analysis

Experiment

Subset Editor

Sample Editor

Analysis

Report

Calculate

Subset: Standards and Unknowns

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|
| A | ● |   |   |   |   |   |   |   |   |    |    |    | ●  | ●  | ●  | ●  | ●  | ●  | ●  | ●  |
| B | ● |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |
| C | ● |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |

Abs Quant results

☒ Positive
 ☒ Negative
 ☒ Uncertain
 ☒ Standard

| Samples                             |       |     |              | Results |               |          |        |
|-------------------------------------|-------|-----|--------------|---------|---------------|----------|--------|
| Include                             | Color | Pos | Name         | Cp      | Concentration | Stand... | Status |
| <input checked="" type="checkbox"/> | ●     |     | A1 Sample 1  | 16.13   | 1.01E5        |          |        |
| <input checked="" type="checkbox"/> | ●     |     | A13 Sample 1 | 16.10   | 1.04E5        |          |        |
| <input checked="" type="checkbox"/> | ●     |     | A14 Sample 1 | 16.13   | 1.02E5        |          |        |
| <input checked="" type="checkbox"/> | ●     |     | A15 Sample 2 |         |               |          |        |
| <input checked="" type="checkbox"/> | ●     |     | A16 Sample 2 | 17.12   | 5.26E4        |          |        |

Export Table


Replicate Statistics

| Samples       | MeanCp | STD Cp | Mean conc | STD conc |
|---------------|--------|--------|-----------|----------|
| A1, A13, A14, | 16.13  | 0.02   | 1.02E5    | 1.44E3   |
| A15, A16, C1, | 17.13  | 0.03   | 5.26E4    | 1.17E3   |
| A17, A18, E1, | 18.27  | 0.02   | 2.47E4    | 2.85E2   |

Apply Template

Notes

Calculate





# AQ Analysis

Experiment

Subset Editor

Sample Editor

**Analysis**

Report




Chart Preferences

Print

Export Chart


Copy to clipboard


# Providing the Standard Curve

- ▶ Use a previously saved standard curve (called an external standard curve). An external standard curve can be loaded into experiments that do not have a standard curve, thus allowing quantitative analysis of those runs. This is especially suitable for applications where the same parameter is analyzed in multiple runs.
- ! *At least one sample (or replicates of this sample) of known concentration must be included in every experiment. This sample should be designated as a standard and should fall within the range of the imported standard curve. The detection format, the analysis mode, and the Color Compensation data (if any) used for the run must be the same as those used for the imported standard curve.*
- ! *For the valid use of the external standard curve, PCR amplification must be highly reproducible and reaction conditions must be constant for all experiments. We recommend running tests to ensure stable PCR efficiency and using replicate samples (especially for low concentrations) to create the standard curve. Also, include a previously quantified sample in each analyzed run, to verify that the calculated values are reproducible.*

## To save a standard curve

On the *Standards* multi-select button, select *Save as external*.



The LightCycler<sup>®</sup> 480 Software automatically navigates to the location *User folder – Special Data – Std Curve* subfolder. To save the standard curve, enter a file name, and click .


You can use the saved external standard curve in other quantification analyses for experiments that have the same experiment parameters as those used to create the standard curve.

# To use an external standard curve


Select *Standard* as the sample type for the standard sample and specify the standard concentration.

For detailed information on the *Sample Editor* see section [Entering Sample Information](#).

On the *Standards* multi-select button, select *Std Curve (External)*.



The *Apply External Standard Curve* dialog opens. Select an appropriate external standard curve object from the list:



# RQ Analysis

77

Experiment

Subset  
EditorSample  
Editor

Analysis

Report

## Create New Analysis

Abs Quant/2nd Derivative Max

Abs Quant/Fit Points

Advanced Relative Quantification

Basic Relative Quantification

## Create new analysis

Analysis Type \* Advanced Relative Quantification

Subset \* All Samples

Program \* amplification

Name \* Advanced Relative Quantification for

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |



## Create new analysis

### Abs Quant Type

☒ Abs Quant/2nd Derivative Max

Sensitivity

☐ High Sensitivity☒ High Confidence☐ Abs Quant/Fit Points

### Subordinate Abs Quant Analysis

☒ Create by Target Name

- Create one analysis for each target name

☐ Create by Filter Combination

- Create one analysis for each filter combination

### Reference Analysis

☒ Create In-Run☐ Select External

### Pairing Rule

☐ One To One☐ All To All☒ All To Mean☐ Mean To All

### Default Standard Curve Settings

When there are no In-Run standards for a target name:

☒ always use efficiency☐ allow external standards with matching target name

# RQ Analysis

Experiment

Subset Editor


Sample Editor

Analysis

Report

Calculate

| Results                             |         |             | Manual Pairing |                    | Target Name |           |            |            |        |
|-------------------------------------|---------|-------------|----------------|--------------------|-------------|-----------|------------|------------|--------|
| Bar                                 |         |             | Target Name    |                    | Target      | Reference | Ratios     |            |        |
| Chart                               | Pairing | Sample Name | Targets        | References         | Mean Cp     | Mean Cp   | Target/Ref | Normalized | Status |
| <input checked="" type="checkbox"/> |         | calibrator  | Target 1       | Reference 1;Refere | 22.43       | 24.14     | 3.266      | 1.000      |        |
| <input checked="" type="checkbox"/> | A1/D1   | sample 1    | Target 1       | Reference 1;Refere | 23.76       | 26.31     | 5.859      | 1.794      |        |
| <input checked="" type="checkbox"/> | A2/D2   | sample 2    | Target 1       | Reference 1;Refere | 21.34       | 24.58     | 9.490      | 2.906      |        |
| <input checked="" type="checkbox"/> | A3/D3   | sample 3    | Target 1       | Reference 1;Refere | 24.30       | 26.18     | 3.675      | 1.125      |        |
| <input checked="" type="checkbox"/> | A4/D4   | sample 4    | Target 1       | Reference 1;Refere | 23.56       | 25.27     | 3.267      | 1.000      |        |
| <input checked="" type="checkbox"/> | A5/D5   | sample 5    | Target 1       | Reference 1;Refere | 25.11       | 23.90     | 0.4320     | 0.1323     |        |
| <input checked="" type="checkbox"/> |         | calibrator  | Target 2       | Reference 1;Refere | 28.49       | 24.14     | 4.91E-2    | 1.000      |        |
| <input checked="" type="checkbox"/> | A7/D1   | sample 1    | Target 2       | Reference 1;Refere | 20.21       | 26.31     | 68.81      | 1402       |        |
| <input checked="" type="checkbox"/> | A8/D2   | sample 2    | Target 2       | Reference 1;Refere | 28.84       | 24.58     | 5.23E-2    | 1.066      |        |
| <input checked="" type="checkbox"/> | A9/D3   | sample 3    | Target 2       | Reference 1;Refere | 33.98       | 26.18     | 4.48E-3    | 9.14E-2    |        |
| <input checked="" type="checkbox"/> | A10/D4  | sample 4    | Target 2       | Reference 1;Refere | 31.32       | 25.27     | 1.51E-2    | 0.3069     |        |
| <input checked="" type="checkbox"/> | A11/D5  | sample 5    | Target 2       | Reference 1;Refere | 33.66       | 23.90     | 1.16E-3    | 2.36E-2    |        |



# RQ Analysis

Experiment

Subset Editor

Sample Editor

Analysis

Report

Settings



## Ratios

☐ Display Target/Reference Ratio

☒ Display Normalized Ratio

☐ Show Ratio Errors

If checked, the ratio shows in the result table and the bar chart.



# Analysis

Experiment

Subset Editor

Sample Editor

Analysis

Report

**Create New Analysis**

- Abs Quant/2nd Derivative Max
- Abs Quant/Fit Points
- Advanced Relative Quantification
- Basic Relative Quantification
- Color Compensation
- Endpoint Genotyping
- Gene Scanning
- Melt Curve Genotyping
- Tm Calling

Create new analysis

**Analysis Type** \* Tm Calling

Subset \* All Samples

**Program** \* Melting Curve

Name \* Tm Calling for All Samples

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| I |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| J |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| K |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| L |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| M |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| N |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| O |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| P |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |

✓
✗



# *T<sub>m</sub> calling/Melting curve*

Experiment



Subset Editor

Sample Editor

Analysis

Report

Calculate



# Report

Experiment

Subset Editor

Sample Editor

Analysis

Report



Report

**Report Settings**

Subset: All Samples

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |

**General** **Detailed**

- ☒ Demo Rel Quant Mono Color
  - ☒ Experiment
  - ☒ Protocol
  - ☐ Samples
  - ☐ Instrument
  - ☐ Revision History
  - ☒ Relative Quantification
    - ☐ Settings
    - ☐ Calibrators
    - ☐ Target Names
    - ☐ Pairings
    - ☐ Results
    - ☐ Results Bar Chart

**Default Settings**

**Apply Template** **Generate**

**LightCycler® 480 Software**

Report

**Demo Rel Quant Mono Color**

Experiment

|               |  |                    |                        |
|---------------|--|--------------------|------------------------|
| Creation Date | 2013/7/24 上午 09:12:10  | Last Modified Date | 2014/3/19 下午 04:44:35  |
| Operator      | Demo   | Owner              | System Admin           |
| Start Time    | 2007/10/15 下午 04:47:45   | End Time           | 2007/10/15 下午 05:58:07 |
| Run State     | Completed  | Software Version   | LCS480 1.4.9.115       |
| Macro         |  | Macro Owner        |                        |
| Macro Status  |  |                    |                        |
| Templates     | Mono Color Hydrolysis Probe - UPL Probe 96-0   | Plate ID           | 00153282               |
| Test ID       |  | Lot ID             |                        |
| Color Comp ID |  |                    |                        |
| Run Notes     | Mono color Relative Quantification experiment. Two target genes and two reference genes are detected with Universal ProbeLibrary probes, all labeled with fluorescein (FAM). |                    |                        |

**Programs**

| Program Name   | Cycles | Analysis Mode  | Target (°C) | Acquisition Mode | Hold (hh:mm:ss) | Ramp Rate (°C/s) | Acquisitions (per °C) | Sec Target (°C) | Step size (°C) | Step Delay (cycles) |
|----------------|--------|----------------|-------------|------------------|-----------------|------------------|-----------------------|-----------------|----------------|---------------------|
| pre-incubation | 1      | None           |             |                  |                 |                  |                       |                 |                |                     |
|                |        |                | 95          | None             | 00:10:00        | 4.40             |                       | 0               | 0              | 0                   |
| amplification  | 45     | Quantification |             |                  |                 |                  |                       |                 |                |                     |
|                |        |                | 95          | None             | 00:00:10        | 4.40             |                       | 0               | 0              | 0                   |
|                |        |                | 60          | Single           | 00:00:30        | 2.20             |                       | 0               | 0              | 0                   |
|                |        |                | 72          | None             | 00:00:01        | 4.40             |                       | 0               | 0              | 0                   |
| cooling        | 1      | None           |             |                  |                 |                  |                       |                 |                |                     |
|                |        |                | 40          | None             | 00:00:30        | 2.20             |                       | 0               | 0              | 0                   |

## Tool Bar



Exit the application



Save changes to selected object



Log into a database



Export the selected object to a file



Show the overview



Close the selected object



Show the navigator



Print the screen



Open Tools



# Create a new detection Format

Tools

- [-] User Access
  - Current Password
  - Users and Groups
  - System Settings
- Report Settings
- Error Log
- [-] Database Information
  - View Logged In Users
  - Update Query Engine
  - Clean-up Database
- Instruments
- Detection Formats

## Detection Formats

| Active                              | Name                        |
|-------------------------------------|-----------------------------|
| <input checked="" type="checkbox"/> | SYBR Green I / HRM Dye      |
| <input checked="" type="checkbox"/> | SimpleProbe                 |
| <input checked="" type="checkbox"/> | Mono Color Hydrolysis Probe |
| <input checked="" type="checkbox"/> | Dual Color Hydrolysis Probe |
| <input checked="" type="checkbox"/> | 3 Color Hydrolysis Probe    |
| <input checked="" type="checkbox"/> | 4 Color Hydrolysis Probe    |
| <input checked="" type="checkbox"/> | Mono Color HybProbe         |
| <input checked="" type="checkbox"/> | Multi Color HybProbe        |

New

Copy

Rename

Delete

## Filter Combination Selection

### Emission

| E | 488 | 510 | 580 | 610 | 640 | 660 |
|---|-----|-----|-----|-----|-----|-----|
| x |     |     |     |     |     |     |
| c |     |     |     |     |     |     |
| i |     |     |     |     |     |     |
| t |     |     |     |     |     |     |
| a |     |     |     |     |     |     |
| t |     |     |     |     |     |     |
| i |     |     |     |     |     |     |
| o |     |     |     |     |     |     |
| n |     |     |     |     |     |     |

Clear

## Selected Filter Combination List

| Excitation Filter | Emission Filter | Name          | Melt Factor | Quant Factor | Max Integration Time (Sec) |
|-------------------|-----------------|---------------|-------------|--------------|----------------------------|
| 465               | 510             | Fluos         | 2           | 1.5          | 2                          |
| 498               | 610             | Red 610       | 1.2         | 5            | 2                          |
| 498               | 640             | Red 640       | 1.2         | 5            | 2                          |
| 498               | 660             | Cy 5 / Cy 5.5 | 1.2         | 5            | 2                          |

Close

Real Time PCR Basic Training

Troubleshooting cases sharing

System Operation Procedures

LC 480 QC Report


Q&A

# LC480 QC Report Template

*Q<sup>+</sup> Qualification Services*  
**Performance Qualification (PQ) Check Report**

*LightCycler® 480 Instrument*

Serial No: 30705  
 Date: 17-Jan, 2022



Roche Diagnostics Ltd.  
 Applied Science Business Area

**Qualification Service**  
 Instrument Verification Run Report

**SYSTEM INFORMATION**

System LightCycler® 480II  
 Instrument Serial No. 30705  
 Instrument SW Version 1.5.1.62  
 Verification Run Started 17-Jan-2022 12:28:44  
 Experiment - File QC-20220117-30705.lxo  
 Reagent Lot Number 57252820  
 Detection - Unit Type II (Filter-Set II)  
 Block - Size 96  
 Comment

**REPORT INFORMATION**

Report Generated With CHECK REPORT TOOL  
 Software Version 8.0.6  
 Parameter Version 1.0

**VERIFICATION RUN RESULT**

Reported Generated By wanga51  
 Report Created 18-Jan-2022  
 Check Result Passed

**Qualification Service**  
 Instrument Verification Run Report

**VERIFICATION RUN REPORT DATA**

| Parameter                        | Lower Limit | Upper Limit | Value  |
|----------------------------------|-------------|-------------|--------|
| <b>Tm Melting Mix (510) TmB2</b> |             |             | Passed |
| Median [°C]                      | 90.3        | 92.7        | 91.42  |
| Range [°C]                       |             | 1.1         | 0.216  |
| <b>Tm Melting Mix (510) dTm2</b> |             |             | Passed |
| Median [°C]                      | -0.35       | 0.25        | -0.074 |
| Range [°C]                       |             | 0.5         | 0.067  |

| Blank Samples       | Excluded | 1 | Passed |
|---------------------|----------|---|--------|
| PCR Mix Samples     |          |   | Passed |
|                     | Excluded | 2 | 0      |
| Melting Mix Samples |          |   | Passed |
|                     | Excluded | 1 | 0      |
| Total Samples       |          |   | Passed |
|                     | Excluded | 3 | 0      |

Document Version 1.0 Page 3/3

**Qualification Service**  
 Instrument Verification Run Report

**VERIFICATION RUN REPORT DATA**

| Parameter                                       | Lower Limit | Upper Limit | Value  |
|---|-------------|-------------|--------|
| <b>Fluorescence 510 offset (Blank Solution)</b> |             |             | Passed |
| Median  | 0.05        | 5.94        | 1.686  |
| CV [%]  |             | 30          | 22.567 |
| <b>Fluorescence 510 background cycle 2-6</b>    |             |             | Passed |
| Mean  | 29.309      | 134.932     | 51.423 |
| CV [%]  |             | 10          | 8.61   |
| <b>Fluorescence 640 signal dynamics</b>         |             |             | Passed |
| Mean  | 0.573       | 1.219       | 0.848  |
| CV [%]  |             | 10          | 8.361  |
| <b>Crossing Point</b>                           |             |             | Passed |
| Mean [Cycles]                                   | 20.9        | 22.9        | 21.882 |
| SD [Cycles]                                     |             | 0.5         | 0.036  |
| <b>Tm PCR Mix (640)</b>                         |             |             | Passed |
| Mean [°C]                                       | 61.1        | 63.5        | 62.409 |
| SD [°C]   |             | 0.3         | 0.053  |
| <b>Tm Melting Mix (510) TmA1</b>                |             |             | Passed |
| Median [°C]                                     | 76.9        | 79.3        | 78.149 |
| Range [°C]                                      |             | 0.9         | 0.121  |
| <b>Tm Melting Mix (510) TmB1</b>                |             |             | Passed |
| Median [°C]                                     | 76.9        | 79.3        | 78.043 |
| Range [°C]                                      |             | 0.9         | 0.333  |
| <b>Tm Melting Mix (510) dTm1</b>                |             |             | Passed |
| Median [°C]                                     | -0.3        | 0.3         | -0.097 |
| Range [°C]                                      |             | 0.5         | 0.255  |
| <b>Tm Melting Mix (510) TmA2</b>                |             |             | Passed |
| Median [°C]                                     | 90.5        | 92.9        | 91.482 |
| Range [°C]                                      |             | 1.1         | 0.168  |

Document Version 1.0 Page 2/3

**CERTIFICATE**  
 FOR INSTRUMENT VERIFICATION RUN

LightCycler® 480II

Verification Run Started: 17-Jan-2022  
 Instrument Serial No.: 30705  
 Certificate Generated By: wanga51  
 Check Report SW Version: 8.0.6  
 Instrument SW Version: 1.5.1.62  
 Reagent Lot: 57252820  
 Experiment - File: QC-20220117-30705.lxo

An instrument Verification Run was carried out with the aforementioned LightCycler® 480II. This certificate confirms that at the time of the verification run this instrument was in accordance with the Roche specifications for the LightCycler® 480II.

Instrument Location: 

Authorized Roche Representative: 

Date: 

Signature: 

Document Version 1.0 1/1

Real Time PCR Basic Training

Troubleshooting cases sharing


System Operation Procedures

LC 480 QC Report

Q&A

# The LightCycler® PRO System

A semi-automated system for research, LDT's and clinical diagnostic workflows



LightCycler®



**Amplify confidence in your results**  
With precise and reproducible data



**Adapt quickly to evolving testing needs**  
With a flexible, semi-open system that supports customisation and standardisation of workflows - based on your lab's unique needs



**Optimise your operational efficiency**  
With convenient, accessible, and reliable data management options

Maximum efficiency with a small footprint



# Roche Digital LightCycler® System

A technological guide to the powerful new addition to our PCR ecosystem



It's time for a leap forward in digital PCR technology. The Roche Digital LightCycler® System is the digital PCR instrument of tomorrow. With a unique combination of 3 nanowell plate configurations, 6 advanced optical channels, and 5x concentrated DNA and RNA master mixes, it has the potential to help your lab to make the leap from publishing research to producing clinically viable assays.



Digital  **LightCycler®**

# Automated Workflows From Sample to Result

**MAGXTRACT 3200**  
 Automated Nucleic Acid purification and PCR Setup System

Beyond Nucleic Acid Purifications

Whirl Mixing


PCR Setup

Open Platform












Magnetic Bead

UV HEPA

Cool Station



Comparison of Bio-Plex Multiplexing vs. ELISA

|  |   |  |  |   |  |
|--|---|--|--|---|--|
| <b>Bio-Plex<br/>Multiplex System</b>                     | <br>1<br>96-well plate     | <br>~3<br>hours       | <br>≤12.5 µl<br>serum or plasma   | <br>50 µl cell<br>culture supernatant      |  |
| <b>What Would it Take to Measure:</b>                    | <br>48 cytokines         | x  | <br>38 samples*                 | =   | <br>1,824 data points |
| <b>Enzyme-Linked<br/>Immunosorbent<br/>Assay (ELISA)</b> | <br>48<br>96-well plates | <br>>106<br>hours** | <br>>1 ml***<br>serum or plasma | <br>>1 ml***<br>cell culture supernatant |  |

提供快速精準的代檢服務  
 Growth factor and Cytokines

\* Samples run in duplicate.  
 \*\* Calculated as 2.2 hours required per plate.  
 \*\*\* Assumes 50 µl of sample used per well.

A large blue hexagon with rounded corners, containing the text "Thank You" in white. The hexagon is set against a background of overlapping light blue and white geometric shapes.

# Thank You

又鑫生物科技有限公司

YU-SHING Bio-Tech Co., Ltd

Sophia Lin

[Sophia.lin@yu-shing.com.tw](mailto:Sophia.lin@yu-shing.com.tw)