



Basic Knowledge of Microscopy

Magi Hung

Classification of Microscopy

NIKON CORPORATION
Instruments Company



Light microscopy

Upright microscopy

Inverted microscopy



Stereomicroscopy



Laser scan microscopy (LSM)



Upright Microscopy

NIKON CORPORATION
Instruments Company



Generally for Slides or sample under shallow medium

Inverted Microscopy


NIKON CORPORATION
Instruments Company



Generally for sample attaching container bottom or slides.
(e.g. cell culture)

Optical train of compound microscope

NIKON CORPORATION
Instruments Company



Nikon CFI60 200/60/25 Spec

NIKON CORPORATION
Instruments Company




200mm : tube lens focal length

60mm : objectives parfocal distance



25mm : diameter of thread

The characteristics of objectives

60x Plan Apochromat Objective



Numerical Aperture (N.A.)




Numerical Aperture = N.A. = $n \cdot \sin \alpha$

α is half the opening angle of the objective.

n is the refractive index of the immersion medium used between the objective and the object.

($n = 1$ for air; $n = 1.51$ for oil or glass)

Resolution



$$d_0 = \frac{1.22\lambda}{N.A._{obj.} + N.A._{Cond}} \quad \text{or more simply } d_0 = \frac{\lambda}{2N.A.}$$

λ = wavelength of light, e.g. 550 nm (green)

● **Resolving power, the limit up to which two small objects are still seen separately.**


Bright Field

Dark Field

Phase Contrast

Differential Interference Contrast

Fluorescence



Hemlock Leaf






- Bright Field is the most universal technique used in light microscope.
- Usually used in samples with colorimetric staining or good contrast.

Dark Field




Fine structures can often not be seen in front of a bright background.

Phase Contrast




Phase Contrast

Living Cells in Brightfield and Phase Contrast




(a)



(b)


Phase Contrast Optical System Alignment



(a)

(c)

(e)



(b)

(d)

(f)




Figure 4



Apodized Phase Contrast

NIKON CORPORATION
Instruments Company




Image with much wider tonal contrast range
See detail like never before

Almost no halo





Previous Phase Image




New Phase (APC) Image

Differential Interference Contrast

NIKON CORPORATION
Instruments Company




HeLa Cell Culture




Heliozoans (*Actinophrys sol*)

Drawback of DIC

Glass vessel only; from cover to bottom.






Vessel cover





When Plastic meets polarized light!

The Principle of Fluorescence




Light path of Fluorescence


NIKON CORPORATION
Instruments Company



Filter Cube




Reduced intensity due to phase plate




- Objective serve as condenser in FL imaging, excitation light output is reduced by phase plate.
- Emission intensity also reduced by phase plate.

Type of Filter Cube




Bandpass emission filter



Longpass emission filter

The Photo-bleaching of Fluorescence

Indian Muntjac deer epidermis fibroblast cells




Nuclei :Hoechst 33258(blue fluorescence)

Mitochondria : MitoTracker Red CMXRos (red fluorescence)


Cytoskeleton :Alexa Fluor 488 (green fluorescence)

Image Digitalization Principle

Sensor Chip




1388



1040



pixel 0,0



pixel 18,22

pixel bitmap



Basic Component of digital image:

- ◆ Pixel (X,Y)
- ◆ Channels (RGB or CYMK)
- ◆ Intensity (12 bit = 0~4096)

Bits does matter!

1 Bit = $2^1 = 2$ steps



4 Bit = $2^4 = 16$ steps



8 Bit = $2^8 = 256$ steps



12 Bit = $2^{12} = 4096$ steps



16 Bit = $2^{16} = 65536$ steps



Note: The human eye can discriminate roughly 200 shades of gray!

Image Format
(Pixel)

Uncompressed File Size
(8 Bit)

256 x 256

64 K

512 x 512

256 K

1024 x 1024

1 MB

2048 x 2048

4 MB

Pixels x bits /8= Bytes/Frame



Nikon 總代理

國祥貿易

turning *vision* into information